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1. Towbin H., Staehlin T. and Gordon J. Electrophoretic transfer of proteins from polyacrylamide gels to nitrocellulose sheets: procedure and some applications. *Proc. Natl. Acad. Sci. U.S.A.* 1976; 76: 4350-4354.
2. Poeszt BJ., Ruscetti FW., Gazdar AF., Bonn PA., Minna J.D. and Gallo RC. Detection and Isolation of type C retrovirus particles from fresh and cultured lymphocytes of a patient with cutaneous T-cell lymphoma. *Proc. Natl. Acad. Sci. U.S.A.* 1980; 77(12): 7415-7419.
3. Kalyanaram VS., Sarngadharan MG., Robert-Guroff M., Miyoshi I., Blayney D., Golde and Gallo R.C. A new subtype of human T-cell leukemia virus (HTLV-II) associated with a T-cell variant of hairy cell leukemia. *Science* 1982; 218: 571-573.
4. William AE., Fang C.T., Slamon D.J. et al. Seroprevalence and epidemiological correlates of HTLV-I infection in U.S. blood donors. *Science* 1988; 240: 643-646.
5. Lee H., Swanson P., Shorty VS., Zack JA., Roseblatt J.D. and Chen I.S.Y. High rate of HTLV-II infection in seropositive IV drug abusers in New Orleans. *Science* 1989; 244: 471-475.

6. Lipka JJ., Bui K., Reyes GR, Moeckli R., Wiktor SZ., Blatner WA., Murphy EL., Hanson CV., Shaw GM., Shinsky JJ. and Fong SKH. Determination of a unique immunodominant epitope of HTLV-I. *Infect Dis* 1990; 162: 353-357

7. Wiktor SZ., Alexandra SS., Shaw GM. et al. Distinguishing between HTLV-I and HTLV-II by Western Blot. *Lancet* 1990; 335: 1533.

8. Samuel KP., Lautenberger JA., Jorcyk CL., Josephs S., Wong Staal F. and Papas TS. Diagnostic potential for human malignancies of bacterially produced HTLV-I envelope protein. *Science* 1984; 226: 1094-1097.

9. Hadlock KG., Goh CJ., Bradshaw PA., Perkins S., Lo J., Habbaz RK., Kaplan J. and Fong SKH. Delineation of an immunodominant and highly HTLV specific epitope within the HTLV-I transmembrane glycoprotein. *Blood* 1995; 68(4): 1392-1399.

10. Varma M., Rudolph D., Knuchel M., Switzer W., Hadlock KG., Velligan M., Chan L., Fong SKH., Lai RB. Enhanced specificity of truncated transmembrane protein for serologic confirmation of HTLV-I and HTLV-II infection by Western Blot assay containing recombinant envelope glycoproteins. *J. Clin. Micro.* 1995; 33(12): 3239-3244.

11. Lillehoj EP, Alexander SS, Dubrule CJ, Wiktor S, Adams R, Thi, A, Manns CC, and Blattner WA. Development and evaluation of a human T-Cell leukemia virus type I serologic confirmatory assay incorporating a recombinant envelope polypeptide. *J. Clin. Microbiol.* 1990; 28: 2653-2658.

12. Lal RB., Brodine SK., Coligan JE., and Roberts CR. Differential antibody responsiveness to p19 gag results in serological discrimination between human T-lymphotropic virus type I and type II. *J. Med. Virol.* 1991; 1: 232-236,

13. Hjelle B., Cyrus S., Swenson S., and Mills R. Serologic distinction between human T-lymphotropic virus (HTLV) type I and HTLV type II. *Transfusion* 1991; 31: 731-736,

14. Madeleine, MM., Wiktor SZ., Goedert JI., Manns A., Levine PH., Biggar RJ., Blattner WA. HTLV-I and HTLV-II worldwide distribution: reanalysis of 4.832 immunoblot results. *Inte. J. Cancer* 1993.; 54(2): 255-260,

15. World Health Organization's Global Programme on AIDS. WHO Global Programme on AIDS Information Update. Virus Information Exchange Newsletter 1990; 7(2): 54-55.

16. Lal, RB., Rudolph DL., Coligan JE., Brodine SK., and Roberts CR. Failure to detect evidence of human T-lymphotropic virus (HTLV) type I and type II in blood donors with isolated gag antibodies to HTLV-I/II. *Blood* 1992; 80: 544-550.

17. Khabbaz, RF., Heneine W., Grindon A., Hartley TM., Shulman G., and Kaplan J. Indeterminate HTLV serologic results in U.S. blood donors: Are they due o HTLV-I or HTLV-II? *J. Acquir. Immune Defic. Syndr.* 1992; 5: 400-404,

18. Lipka, JJ., Young KK., Kwok SY., Reyes GR., Sninsky JJ, and Fong SK. Significance of human T-lymphotropic virus type I indeterminate serological findings among healthy individuals. *Vox Sang.* 1991; 61: 171-176,

19. Zrein M., Louwagie J., Boeykens H., Govers L., Hendicks G., Bosman F., Sablon E., Demarquilly C., Boniface M., and Saman E. Assessment of a new immunoassay for serological confirmation and discrimination of human T-cell lymphotropic virus infections. *Clin. Diag. Lab. Immunol.* 1998; 5: 45-49.

20. Witt DJ., Kuramoto K., Kemper M., and Holland P. Utility of prospective study of donors deferred as HTLV indeterminate. *Vox Sang.* 2000; 78:130-131,

21. Hayes C.G., Burans J.P., and Oberst R.B. Antibodies to Human T Lymphotropic virus Type I in a population from the Philippines: Evidence for cross-reactivity with *Plasmodium falciparum*. The J. Infect. Dis. 1990; 163: 257-262,


22. Gallo D., Diggs JL., and Hanson CV. Evaluation of two commercial Human T-Cell Lymphotropic Virus Western blot (Immunoblot) kits with problem specimens. *J. Clin. Microbiol.* 1994; 32: 2046-2049.

23. Garin B., Gosselin S., de The G., and Gessain A. HTLV-II infection in a high viral endemic area of Zaire, Central Africa: Comparative evaluation of serology, PCR, and significance of indeterminate Western blot pattern. *J. Med. Virol.* 1994; 44: 104-109.

24. Fujiyama C., Fujiyoshi T., Matsumoto D., Yashiki S., Tamashiro H., and Sonoda S. Re-evaluation of anti-HTLV-I Western blot assay using HTLV-I and HTLV-II serum panels. Clin. & Diag. Virol. 1995; 4: 149-161,

25. Rouet F, Meertens L, Courouble G, Herrmann-Storck C, Pabinger U, Chancellor B, Abid A, Strobel M, Maulcuere P, and Gessain A. Serological, epidemiological, and molecular differences between human T-cell lymphotropic virus type I - seropositive healthy carriers and persons with HTLV-I gag indeterminate Western blot patterns from the Caribbean. *J. Clin. Microbiol.* 2001; 39: 1247-1253
26. Cesaire R., Bera O., Maier H., Lezin A., Martial J., Ouka M., Kerob-Bauchet B., Ould Amar AK., and Vernant JC. Seroindeterminate patterns and seroconversions to human T-lymphotropic virus type I positivity in blood donors from Martinique, French West Indies. *Transfusion* 1999; 39: 1145-1149.

27. Soldan SS., Graf MD., Waziri A., Flerlage AN., Robinson SM., Kawaninshi T., Leist TP., Lehky T.J., Levin MC., and Jacobson S. HTLV-I/II seroindeterminate western blot reactivity in a cohort of patients with neurological disease. *J. Infect. Dis.* 1999; 180: 685-694.

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